

TECHNICAL DATA SHEET

SteriCleen[®] Scale Remover 5L

Acid-based Medical Instrument & Equipment Cleaner Australian Medical Device Class 1 ARTG No. 332455

DESCRIPTION:

SteriCleen[®] Scale Remover is a highly active acid medical device detergent specifically designed for the removal of oxides, salts and other inorganic residues via both machine and manual cleaning operations. SteriCleen[®] Scale Remover is a concentrated formula which is low foaming, non-staining and effective in both hard and soft water.

BENEFITS:

- · Concentrated acid solution for fast descaling
- · Formulated to be safe for use on most instruments and devices
- Biodegradable surfactant system with recyclable packaging
- Phosphate-free
- Non-toxic, non-abrasive, low foaming, free rinsing.
- Contains corrosion inhibitors for the prevention of corrosion of medical instruments during the cleaning process.
- Can be used in Ultrasonic & Manual cleaning of instruments prior to disinfecting or sterilisation procedures.
- Complies with AS/NZS4187 and AS/NZS4815

INTENDED USAGE:

- Manual cleaning of Medical and Dental instruments prior to the disinfection or sterilisation process.
- Use as a dedicated cleaner in conjunction with Ultrasonic cleaning machines.
- Manual cleaning of large and small equipment.

DIRECTIONS FOR USE:

- 1. Manufacturer's instructions should be strictly followed as advised by the instrument provider.
- 2. Prepare a fresh ready to use solution daily.
- 3. Dilution prior to use:
 - Heavy soiling: Dilute 20ml per Litre of water
 - Medium to light soiling: Dilute 10mL per litre of water

Discard diluted solutions within 24 hours of preparation

Note: Instrument cleaners are not disinfectants, and all instruments should be cleaned, rinsed and drained before immersion in an instrument grade DISINFECTANT, or prior to autoclaving.

CODE: 9224/42







TECHNICAL DATA SHEET

CHEMICAL COMPOSITION AND PROPERTIES:

Colour concentrate:
Colour dilute solution:
Odour (Concentrate):
pH concentrate:
pH dilution:
Solubility:
Rinsibility:
Specific gravity:
Corrosion inhibitor:
Shelf life:

Green Colourless No odour 2.5 + 0.5 2.5 ± 1.0 Completely soluble in water Residue free 1.02 ± 0.05 g/ml Proprietary blend 3 years

STORAGE CONDITIONS: Store below 25°C. Do not expose to direct sunlight.

SAFETY DATA: Information is listed in the Material Safety Data Sheet.

PACK SIZE: 5 Litres

WARRANTY – All statements, information and data presented herein are believed to be accurate and reliable but are not to be taken as a guarantee, expressed or implied, for which seller assumes legal responsibility and they are offered solely for your consideration, investigation and verification. Statements or suggestions concerning possible use of this product are made without representation or warranty that any such use is free of patent infringement and are not recommendations to infringe on any patent. Created 30 November 2021 Date Printed 1/08/2022 11:07 AM