

## TECHNICAL DATA SHEET

# SteriCleen Multi-Enzymatic Neutral Detergent 5L

CODE: 9227/42

Neutral Medical Instrument & Equipment Cleaner  
Australian Medical Device Class 1  
ARTG No. 332455

### DESCRIPTION:

SteriCleen Multi-Enzymatic Neutral Detergent is a highly active neutral multi-enzyme medical device detergent specifically formulated for the removal of hospital and surgical soils such as blood, body tissue, polysaccharides, lipids and other protein-based soils via both machine and manual cleaning operations. SteriCleen is a concentrated formula which is low-foaming, non-staining and effective in both hard and soft water.

### BENEFITS:

- Concentrated blend with neutral pH– safe for use on most instruments and devices
- Biodegradable surfactant system with recyclable packaging
- Acts quickly on blood and other protein-based soils
- Non-toxic, non-abrasive, low foaming, free rinsing
- Contains corrosion inhibitors for the prevention of corrosion of medical instruments during the cleaning process.
- Can be used in Ultrasonic & Manual cleaning of instruments prior to disinfecting or sterilisation procedures.
- Complies with AS/NZS4187 and AS/NZS4815

### INTENDED USAGE AND DIRECTIONS TO USE:

- Use as a dedicated cleaner in conjunction with Ultrasonic cleaning machines.
- Manual cleaning of Medical and Dental instruments prior to the disinfection or sterilisation process.
- Manual cleaning of theatre floors, tables, benches as well as laboratory glassware.
- Manual cleaning of large and small equipment.

### DIRECTIONS FOR USE:

1. Manufacturer's instructions should be strictly followed as advised by the instrument provider.
2. Prepare a fresh ready to use solution daily.
3. Dilution instructions:

Soil and contamination level	Dilution Rate	Recommended soak time
HEAVY	1:100 (10mL per litre)	7 – 9 minutes
MEDIUM	1:150 (7.5mL per litre)	7 minutes
LIGHT	1:200 (5mL per litre)	5 minutes



## TECHNICAL DATA SHEET

Use warm water to dilute to achieve optimum performance (30° C - 40° C).

Do not allow biological contaminants to dry. If Contaminants have dried, soak at least for 20 mins at a dilution of 1:50.

Rinse with water after removal from soak solution and prior to sterilisation / disinfection.

4. Discard solution after use.

Note: Instrument cleaners are not disinfectants, and all instruments should be cleaned, rinsed, and drained before immersion in an instrument grade DISINFECTANT, or prior to autoclaving.

### CHEMICAL COMPOSITION AND PROPERTIES

Colour concentrate:	Light Blue
Colour dilute solution:	Yellow
Odour (Concentrate):	no odour
pH concentrate:	7.0-9.0
pH dilution:	7.5 ± 0.5
Solubility:	Completely soluble in water
Rinsibility:	Residue free
Specific gravity:	1.02 ± 0.05g/ml
Corrosion inhibitor:	Proprietary blend
Shelf life:	1 year

**STORAGE CONDITIONS:** Store below 25°C. Do not expose to direct sunlight.

**SAFETY DATA:** Information is listed in the Material Safety Data Sheet.

**PACK SIZE:** 5 Litres

## TECHNICAL DATA SHEET

**WARRANTY** – All statements, information and data presented herein are believed to be accurate and reliable but are not to be taken as a guarantee, expressed or implied, for which seller assumes legal responsibility and they are offered solely for your consideration, investigation, and verification. Statements or suggestions concerning possible use of this product are made without representation or warranty that any such use is free of patent infringement and are not recommendations to infringe on any patent.

Created 30

November 2021 Date Printed 23/05/2022 12:15 AM